

Breast Examination Trainer - Standard, Advanced & Basic

40200 Standard


40201 Advanced

40202 Basic

Simulated Patient



Breasts can be worn by a simulated patient (male or female)



Pathologies are supplied and fit securely into the holes in the Pathology Supports & Back Plates

Lymph Node Pads

Torso

Our Breast Examination Trainer provides a highly realistic learning platform for acquiring the skills required to perform Clinical Breast Examination (CBE). Featuring 6 readily interchangeable and multi-positional pathologies, providing healthcare professionals with the tools to identify various complications and pathologies, including carcinomas, cysts, fibrocystic disease and fibroadenoma. Both Simulated Patient and bench top training can be used for any undergraduate programs

running OSCEs or healthcare professionals promoting best practice techniques to trainees and patients.

This product was developed in close collaboration with leading clinicians:

- Mr Peter Donnelly, Torbay Hospital, UK
- Clinical Skills Resources Centre, University of Liverpool, UK
- Clinical Practice Centre, St James's University Hospital, Leeds, UK

Product Qualities

Anatomy

- Realistic soft tissue breast anatomy
- Pathologies supplied: carcinomas: 2cm, 3cm, 5cm, cyst, fibrocystic disease, fibroadenoma

Realism

- Soft tissue breasts look and feel realistic
- Clavicular and axilla pads for accurate lymph node placement

Key Features

- Pathologies can be placed in various predetermined location points and are easily changeable
- Dual purpose product: bench top and Simulated Patient
- Hard torso available for bench top use

Versatility

- Comfortable to wear for long periods of time during OSCEs and assessments

Safety

- Latex free

Cleaning & Care

- Skin surface is washable using soap and water

Skills

- Clinical Breast Examination (CBE)
- Self Breast Examination (SBE)
- Identification of anatomical landmarks
- Identification of lymph nodes (axillary, supra and infraclavicular)
- Location and diagnosis of pathologies
- Professional-to-patient communication

	Adv	Std	Basic
• Wearable Examination Breasts	✓	✓	✓
• 40204 Pathologies	✓	✓	✓
• 40205 Breast Examination Inserts	✓	✓	✓
• 40206 Breast Pathologies Supports	✓	✓	✓
• 40207 Breast Back Plates	✓	✓	✓
• 40208 Breast Examination Torso	✓	✓	✗
• 40203 Lymph Node Pads + Lymph Nodes	✓	✗	✗

Optional Extras

- **40209** Breast Examination Carry Case

Patents granted

• US6336812 • EP099027 • DE69804793 • CA 2293585
• ZL98806331 • HK1023832 • 755575 Australia