

Venepuncture Arms

- Standard
- Advanced
- Bag Stand & Advanced



Learn more at

limbsandthings.com

**Learning Has Never
Been More Lifelike**

Product Range

Standard Venepuncture Arm:



The Standard Venepuncture Arm is ideal for a clinical skills training environment where large numbers of trainees are practising this skill on a regular basis, and for student self-directed learning, due to ease of set up and maintenance. It comes with a Bag & Stand Mock Blood Supply and has the following accessible veins: cephalic, basilic and 2 dorsal metacarpals.

Advanced Venepuncture Arm:



The Advanced Venepuncture Arm has a greater number of accessible veins: cephalic, basilic, median cubital and 4 dorsal metacarpals. It comes with an adjustable pressurized fluid system. It addresses the fundamental needs of clinical skills tutors, trainees and technicians alike.

Bag Stand & Advanced Venepuncture Arm:



This product combines the vein structure of the Advanced Venepuncture Arm with the simplicity of the Bag & Stand Mock Blood Supply.

Skills Gained

- Venepuncture
- IV Cannulation
- IV Infusion

Product Qualities



Standard Venepuncture Arm:

Part No:

- 00330 Light
- 00331 Brown
- 00332 Black



You'll Receive:

- 00333 1 SV Arm Vein
- 00334 1 SV Vein Module
- 00336 1 Bag & Stand Venepuncture Arm Mock Blood Supply
- 00020 1 Concentrated Venous Blood - Starter Pack
- 00335 1 Venepuncture Arm Infusion Tube
- 00302 1 Venepuncture Arm Rest
- 1 Venepuncture Arm Skin
- 3 Sachets Water-based Lubricant
- 1 Arm Shell

Packed weight (kg) & dimensions (mm) 6.1 kg 1010 l 190 w 290 h

Advanced Venepuncture Arm:

Part No:

- 00290 Light
- 00298 Brown
- 00296 Black



You'll Receive:

- 00292 1 AV Arm Vein
- 00295 1 AV Vein Module
- 00020 1 Concentrated Venous Blood - Starter Pack
- 00335 1 Venepuncture Arm Infusion Tube
- 00302 1 Venepuncture Arm Rest
- 00303 1 Refill Bottle
- 00337 1 Pressurised Venepuncture Arm Mock Blood Supply
- 1 Venepuncture Arm Skin
- 3 Sachets Water-based Lubricant
- 1 Arm Shell

Packed weight (kg) & dimensions (mm) 4.25 kg 1010 l 190 w 290 h

Bag Stand & Advanced Venepuncture Arm:

Part No:

- 00440 Light
- 00441 Brown
- 00442 Black



You'll Receive:

- 00292 1 AV Arm Vein
- 00295 1 AV Vein Module
- 00336 1 Bag & Stand Venepuncture Arm Mock Blood Supply
- 00020 1 Concentrated Venous Blood - Starter Pack
- 00335 1 Venepuncture Arm Infusion Tube
- 00302 1 Venepuncture Arm Rest
- 1 Venepuncture Arm Skin
- 3 Sachets Water-based Lubricant
- 1 Arm Shell

Packed weight (kg) & dimensions (mm) 6.1 kg 1010 l 190 w 290 h

+ Optional Extras & Consumables: Available for all Venepuncture Arms

- 00291 Arm Skin - Light
- 00299 Arm Skin - Brown
- 00297 Arm Skin - Black
- 00300 Venepuncture Arm Harness
- 00301 Venepuncture Arm Carry Case



Learn more at:

limbsandthings.com

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501

E: UK Customers - gbsales@limbsandthings.com

E: International Customers - sales@limbsandthings.com